

WorkstationJ Keyboard Shortcuts (September 2010)

Keyboard Shortcuts for the **MultiModality Viewer Three Registered CT** Plugin:

1. /: go to next case
2. Left Arrow (←) : go to previous image in current case
3. Right Arrow (→) : go to next image in current case
4. Space Bar : go to next image in current case
5. Enter: go to next image in current case
6. Asterisk (*): show/hide crosshairs

Keyboard Shortcuts for the **Single-CT Chest Viewer** Plugin:

1. /: go to next case
2. Plus (+): increase magnification
3. Minus (-) decrease magnification